

PKC gamma Rabbit mAb


ABZ3121

Product Overview


Name	PKC gamma Rabbit mAb
Accession(Primary)	P05129
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P05129)

Gene Name	PRKCG
About Gene Symbol	PRKCG (Protein Kinase C Gamma) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PKC gamma in rat brain lysates using PKC gamma antibody.



Western blot analysis of PKC in Hela lysates using PKC antibody.