

PMCA1 Rabbit mAb


ABZ3122

Product Overview


Name	PMCA1 Rabbit mAb
Accession(Primary)	P20020
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 135 kDa; Observed MW: 130-250 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; IP-1:20-1:50

Target information (P20020)

Gene Name	ATP2B1
About Gene Symbol	



Western blot analysis of PMCA1 in mouse brain, rat lung lysates using PMCA1 antibody.



Western blot analysis of PMCA1 in U87-MG, Human placenta lysates using PMCA1 antibody.