


PON1 Rabbit mAb

ABZ3125


Product Overview

Name	PON1 Rabbit mAb
Accession(Primary)	P27169
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20

Target information (P27169)

Gene Name PON1**About Gene Symbol**

Western blot analysis of PON1 in mouse liver lysates using PON1 antibody.



Western blot analysis of PON1 in K562 lysates
using PON1 antibody.