

PPM1G Rabbit mAb


ABZ3126

Product Overview

Name	PPM1G Rabbit mAb
Accession(Primary)	O15355
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:200

Target information (O15355)

Gene Name	PPM1G
About Gene Symbol	PPM1G (Protein Phosphatase, Mg ²⁺ /Mn ²⁺ Dependent 1G) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PPM1G in K562, rat Brain, C6, 3T3, HeLa lysates using PPM1G antibody.