

PRC1 Rabbit mAb


ABZ3127

Product Overview


Name	PRC1 Rabbit mAb
Accession(Primary)	O43663
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:20-1:50; FC-1:20-1:50; IHC-1:50-1:100

Target information (O43663)

Gene Name	PRC1
About Gene Symbol	PRC1 (Protein Regulator Of Cytokinesis 1) is a cytokine and immune signaling protein involved in JAK-STAT/cytokine signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PRC1 in C6, Jurkat lysates using PRC1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PRC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.