

Prealbumin Rabbit mAb


ABZ3128

Product Overview

Name	Prealbumin Rabbit mAb
Accession(Primary)	P02766
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:10-1:100

Target information (P02766)

Gene Name	TTR
About Gene Symbol	



Western blot analysis of Prealbumin in Human tonsil, Human placenta lysates using Prealbumin antibody.