

Prion Protein PrP Rabbit mAb


ABZ3129

Product Overview


Name	Prion Protein PrP Rabbit mAb
Accession(Primary)	P04156
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:20-1:50

Target information (P04156)


Gene Name	PRNP
About Gene Symbol	



Western blot analysis of Prion Protein in U251 lysates using Prion Protein antibody.



Western blot analysis of Prion protein in rat Brain and mouse Brain lysates using Prion Protein antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Prion Protein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.