

PRKACG Rabbit mAb

ABZ3130


Product Overview

Name	PRKACG Rabbit mAb
Accession(Primary)	P22612
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000


Target information (P22612)

Gene Name PRKACG

About Gene Symbol



Western blot analysis of PRKACG in mouse testis, mouse lung lysates using PKA C gamma antibody.



Western blot analysis of PRKACG in MCF-7, C6, HeLa lysates using PRKACG antibody.