

PRMT3 Rabbit mAb


ABZ3131

Product Overview


Name	PRMT3 Rabbit mAb
Accession(Primary)	O60678
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:20-1:50; IP-1:20-1:50

Target information (O60678)

Gene Name	PRMT3
About Gene Symbol	PRMT3 (Protein Arginine Methyltransferase 3) is a methyltransferase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PRMT3 in rat Brain, C6, 3T3, Hela lysates using PRMT3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PRMT3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.