

PRKAR1A Rabbit mAb

ABZ3133


Product Overview

Name	PRKAR1A Rabbit mAb
Accession(Primary)	P10644
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P10644)

Gene Name PRKAR1A

About Gene Symbol



Western blot analysis of Protein Kinase A regulatory subunit I alpha/PRKAR1A in C6, rat Brain, Jurkat lysates using PKA RI alpha antibody.

