

Prothrombin Rabbit mAb

ABZ3134


Product Overview

Name	Prothrombin Rabbit mAb
Accession(Primary)	P19221
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-:1:200-1:2000

Target information (P19221)

Gene Name F2

About Gene Symbol



Western blot analysis of Prothrombin in mouse liver, rat liver lysates using Prothrombin antibody.

