

PRPF4 Rabbit mAb


ABZ3135

Product Overview

Name	PRPF4 Rabbit mAb
Accession(Primary)	O43172
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:500-1:1000; FC-1:20-1:50; IP-1:20-1:50

Target information (O43172)

Gene Name	PRPF4
About Gene Symbol	PRPF4 (Pre-mRNA Processing Factor 4) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PRPF4 in K562, 3T3, HeLa lysates using PRPF4 antibody.