


PRSS8 Rabbit mAb

ABZ3136

Product Overview

Name	PRSS8 Rabbit mAb
Accession(Primary)	Q16651
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q16651)

Gene Name PRSS8**About Gene Symbol**

Western blot analysis of PRSS8 in MCF-7 lysates
using PRSS8 antibody.

