

PTN Rabbit mAb

ABZ3137


Product Overview

Name	PTN Rabbit mAb
Accession(Primary)	P21246
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:500-1:1000


Target information (P21246)

Gene Name

PTN

About Gene Symbol

Western blot analysis of PTN in U87-MG lysates using PTN antibody.



Western blot analysis of PTN in C6, rat Brain
lysates using PTN antibody.