

PTP1B Rabbit mAb


ABZ3138

Product Overview

Name	PTP1B Rabbit mAb
Accession(Primary)	P18031
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:20-1:50

Target information (P18031)

Gene Name	PTPN1
About Gene Symbol	



Western blot analysis of PTP1B in 293T, HCT116
lysates using PTP1B antibody.