


Rad51D Rabbit mAb

ABZ3140

Product Overview

Name	Rad51D Rabbit mAb
Accession(Primary)	O75771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; IP-1:20-1:50

Target information (O75771)

Gene Name RAD51D**About Gene Symbol**

Western blot analysis of Rad51D in K562, C6, 3T3, HeLa lysates using Rad51D antibody.

