


RAMP1 Rabbit mAb

ABZ3142

Product Overview

Name	RAMP1 Rabbit mAb
Accession(Primary)	O60894
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O60894)

Gene Name RAMP1**About Gene Symbol**

Western blot analysis of RAMP1 in Raji, COS7, C6, 3T3, HeLa lysates using RAMP1 antibody.

