

RAVER2 Rabbit mAb


ABZ3143

Product Overview



Name	RAVER2 Rabbit mAb
Accession(Primary)	Q9HCJ3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; IP-1:10-1:100

Target information (Q9HCJ3)


Gene Name	RAVER2
About Gene Symbol	RAVER2 (Ribonucleoprotein, PTB Binding 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of RAVER2 (green) in K562 using RAVER2 antibody, and DAPI (blue).



Western blot analysis of RAVER2 in Jurkat, HeLa lysates using RAVER2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using RAVR2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.