

## Rb Rabbit mAb


ABZ3144

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Rb Rabbit mAb                                |
| <b>Accession(Primary)</b>   | P06400                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 106 kDa; Observed MW: 110 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                              |

### Target information (P06400)

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | RB1                                                                                                                                                                                                                                                           |
| <b>About Gene Symbol</b> | RB1 (RB Transcriptional Corepressor 1) is a protein phosphatase involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Rb in 3T3 lysates using Rb antibody.