

REEP2 Rabbit mAb


ABZ3145

Product Overview

Name	REEP2 Rabbit mAb
Accession(Primary)	Q9BRK0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9BRK0)

Gene Name	REEP2
About Gene Symbol	REEP2 (Receptor Accessory Protein 2) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of REEP2 in Raoms, Jurkat
lysates using REEP2 antibody.