

Relaxin 2/RLN2 Rabbit mAb

ABZ3146


Product Overview


Name	Relaxin 2/RLN2 Rabbit mAb
Accession(Primary)	P04090
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (P04090)

Gene Name RLN2

About Gene Symbol





Immunohistochemistry analysis of paraffin-embedded Human tonsil using Relaxin 2/RLN2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.