

Retinoic Acid Receptor beta Rabbit mAb


ABZ3147

Product Overview

Name	Retinoic Acid Receptor beta Rabbit mAb
Accession(Primary)	P10826
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 50 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P10826)

Gene Name	RARB
About Gene Symbol	RARB (Retinoic Acid Receptor Beta) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Retinoic Acid Receptor
beta in U87-MG lysates using Retinoic Acid
Receptor beta antibody.