

RPL10A Rabbit mAb


ABZ3151

Product Overview

Name	RPL10A Rabbit mAb
Accession(Primary)	P62906
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:20-1:50

Target information (P62906)

Gene Name	RPL10A
About Gene Symbol	



Western blot analysis of RPL10A in K562, C6, 3T3, HeLa lysates using RPL10A antibody.