

RRM1 Rabbit mAb


ABZ3153

Product Overview


Name	RRM1 Rabbit mAb
Accession(Primary)	P23921
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 90 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (P23921)


Gene Name	RRM1
About Gene Symbol	RRM1 (Ribonucleotide Reductase Catalytic Subunit M1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms.



Western blot analysis of RRM1 in U2OS, HepG2, MCF-7 lysates using RRM1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using RRM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of RRM1 (green) in Hela using RRM1 antibody, and DAPI (blue)