

RRP4 Rabbit mAb


ABZ3154

Product Overview


Name	RRP4 Rabbit mAb
Accession(Primary)	Q13868
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; FC-1:100-1:200

Target information (Q13868)

Gene Name	EXOSC2
About Gene Symbol	



Western blot analysis of RRP4 in rat spleen lysates using RRP4 antibody.



Western blot analysis of RRP4 in 3T3, HeLa lysates using RRP4 antibody.