

SEC6 Rabbit mAb


ABZ3155

Product Overview


Name	SEC6 Rabbit mAb
Accession(Primary)	O60645
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 86 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; IHC-1:50-1:100

Target information (O60645)

Gene Name	EXOC3
About Gene Symbol	



Western blot analysis of rSec6 in rat Brain, 3T3, HeLa lysates using Exocyst Complex Component 3 antibody.



Immunocytochemistry analysis of rSec6(green) in HeLa using rSec6 antibody, and DAPI(blue)