

S6K1 Rabbit mAb


ABZ3156

Product Overview


Name	S6K1 Rabbit mAb
Accession(Primary)	P23443
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:100-1:2000; IP-1:50-1:100

Target information (P23443)


Gene Name	RPS6KB1
About Gene Symbol	



Western blot analysis of S6K1 in Raw264.7, 3T3 lysates using S6K1 antibody.



Western blot analysis of S6K1 in rat Brain lysates using S6K1 antibody.



Immunocytochemistry analysis of S6K1(green) in
Hela using S6K1 antibody, and DAPI(blue)