

SAMD9 Rabbit mAb


ABZ3157

Product Overview

Name	SAMD9 Rabbit mAb
Accession(Primary)	Q5K651
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 184 kDa; Observed MW: 184 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (Q5K651)

Gene Name	SAMD9
About Gene Symbol	



Western blot analysis of SAMD9 in A549 lysates
using SAMD9 antibody.