

SCD1 Rabbit mAb


ABZ3158

Product Overview


Name	SCD1 Rabbit mAb
Accession(Primary)	O00767
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:400; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (O00767)


Gene Name	SCD
About Gene Symbol	



Western blot analysis of SCD1 in HeLa lysates using SCD1 antibody.



Immunocytochemistry analysis of SCD1(green) in HeLa using SCD1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SCD1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.