

PLA2G2A Rabbit mAb

ABZ3159


Product Overview

Name	PLA2G2A Rabbit mAb
Accession(Primary)	P14555
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P14555)

Gene Name PLA2G2A

About Gene Symbol



Western blot analysis of Secretory Phospho-lipase A2 in K562, C6, 3T3, HeLa lysates using Secretory Phospho-lipase A2 antibody.

