

SGSM1 Rabbit mAb


ABZ3162

Product Overview


Name	SGSM1 Rabbit mAb
Accession(Primary)	Q2NKQ1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 130 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (Q2NKQ1)

Gene Name	SGSM1
About Gene Symbol	



Western blot analysis of SGSM1 in K562, HeLa lysates using SGSM1 antibody.



Immunocytochemistry analysis of SGSM1(green) in HeLa using SGSM1 antibody, and DAPI(blue)