


SH2D1A/SAP Rabbit mAb

ABZ3163


Product Overview

Name	SH2D1A/SAP Rabbit mAb
Accession(Primary)	O60880
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O60880)

Gene Name SH2D1A**About Gene Symbol**

Western blot analysis of SH2D1A/SAP in 3T3
lysates using SH2D1A antibody.



Western blot analysis of SH2D1A/SAP in MCF-7, C6 lysates using SH2D1A/SAP antibody.