

SH3PX1/SNX9 Rabbit mAb


ABZ3164

Product Overview


Name	SH3PX1/SNX9 Rabbit mAb
Accession(Primary)	Q9Y5X1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 67,78 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:10-1:100; IP-1:20-1:50

Target information (Q9Y5X1)


Gene Name	SNX9
About Gene Symbol	SNX9 (Sorting Nexin 9) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of SNX9 in mouse intestines lysates using SNX9 antibody.



Western blot analysis of SH3PX1 in K562, HeLa lysates using SH3PX1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using SNX9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of SNX9 (green) in Hela using SNX9 antibody, and DAPI (blue).