


SODD Rabbit mAb

ABZ3165

Product Overview

Name	SODD Rabbit mAb
Accession(Primary)	O95429
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 50 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (O95429)

Gene Name BAG4**About Gene Symbol**

Western blot analysis of Silencer of Death Domain
in MOLT4 lysates using SODD antibody.

