

SLBP Rabbit mAb


ABZ3166

Product Overview

Name	SLBP Rabbit mAb
Accession(Primary)	Q14493
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:50; IP-1:20

Target information (Q14493)

Gene Name	SLBP
About Gene Symbol	SLBP (Stem-Loop Histone mRNA Binding Protein) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics and translational regulation. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SLBP in 293T, Jurkat, Ramos lysates using SLBP antibody.