

SLC25A12 Rabbit mAb


ABZ3167

Product Overview

Name	SLC25A12 Rabbit mAb
Accession(Primary)	O75746
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 75 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (O75746)

Gene Name	SLC25A12
About Gene Symbol	SLC25A12 (Solute Carrier Family 25 Member 12) is a membrane transporter involved in membrane transport. It is an important research target in pharmacology, drug resistance, and membrane biology and is widely studied to elucidate cellular mechanisms and id



Western blot analysis of SLC25A12 in rat Brain, C6, 3T3 lysates using SLC25A12 antibody.