

SNAP29 Rabbit mAb


ABZ3168

Product Overview



Name	SNAP29 Rabbit mAb
Accession(Primary)	O95721
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:10-1:100; IP-1:20-1:50

Target information (O95721)

Gene Name	SNAP29
About Gene Symbol	SNAP29 (Synaptosome Associated Protein 29) is a membrane receptor involved in vesicular trafficking. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of SNAP29 (green) in K562 using SNAP29 antibody, and DAPI (blue).



Western blot analysis of SNAP29 in rat Brain, 3T3, HeLa lysates using SNAP29 antibody.