

SMARCA5 Rabbit mAb

ABZ3169


Product Overview

Name	SMARCA5 Rabbit mAb
Accession(Primary)	O60264
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 122 kDa; Observed MW: 122 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (O60264)


Gene Name SMARCA5

About Gene Symbol



Immunocytochemistry analysis of SMARCA5 (green) in Jurkat using SMARCA5 antibody, and DAPI(blue).





Western blot analysis of SNF2H in K562, C6, 3T3, HeLa lysates using SNF2H antibody.