

SOCS2 Rabbit mAb


ABZ3172

Product Overview

Name	SOCS2 Rabbit mAb
Accession(Primary)	O14508
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:10-1:100

Target information (O14508)

Gene Name	SOCS2
About Gene Symbol	



Western blot analysis of SOCS2 in K562 lysates
using SOCS2 antibody.