


Somatostatin Receptor 3 Rabbit mAb

ABZ3173

Product Overview

Name	Somatostatin Receptor 3 Rabbit mAb
Accession(Primary)	P32745
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:1000

Target information (P32745)

Gene Name SSTR3**About Gene Symbol**

Western blot analysis of Somatostatin Receptor 3/SSTR3 in Jurkat lysates using Somatostatin Receptor 3 antibody.

