

## SP1 Rabbit mAb


ABZ3174

### Product Overview


|                             |                                                |
|-----------------------------|------------------------------------------------|
| <b>Name</b>                 | SP1 Rabbit mAb                                 |
| <b>Accession(Primary)</b>   | P08047                                         |
| <b>Host Species</b>         | Rabbit                                         |
| <b>Clonality</b>            | Monoclonal                                     |
| <b>Isotype</b>              | IgG                                            |
| <b>Species Reactivity</b>   | Human, Mouse                                   |
| <b>Calculated MW</b>        | Calculated MW: 81 kDa; Observed MW: 90 kDa     |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200 |

### Target information (P08047)


|                          |                                                                                                                                                                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | SP1                                                                                                                                                                                                                                                                                      |
| <b>About Gene Symbol</b> | SP1 (Sp1 Transcription Factor) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Immunocytochemistry analysis of SP1 (green) in Hela using SP1 antibody, and DAPI (blue).



Western blot analysis of SP1 in K562, 3T3, Hela lysates using SP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using SP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.