

SP1 Rabbit mAb


ABZ3175

Product Overview


Name	SP1 Rabbit mAb
Accession(Primary)	P08047
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 81 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1; 50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P08047)

Gene Name	SP1
About Gene Symbol	SP1 (Sp1 Transcription Factor) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using SP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of SP1 (green) in Hela using SP1 antibody, and DAPI (blue)

