


SPOP Rabbit mAb

ABZ3177

Product Overview

Name	SPOP Rabbit mAb
Accession(Primary)	O43791
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O43791)

Gene Name SPOP**About Gene Symbol**

Western blot analysis of SPOP in MCF-7 lysates using SPOP antibody.

