

Pulmonary Surfactant Associated Protein D Rabbit mAb


ABZ3179

Product Overview

Name	Pulmonary Surfactant Associated Protein D Rabbit mAb
Accession(Primary)	P50404
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:1000-1:5000

Target information (P50404)

Gene Name	Sftpd
About Gene Symbol	



Western blot analysis of Surfactant protein D in
rat lung lysates using Pulmonary Surfactant
Associated Protein D antibody.