

# TROP2 Rabbit mAb


ABZ3183

## Product Overview

|                             |                                                                          |
|-----------------------------|--------------------------------------------------------------------------|
| <b>Name</b>                 | TROP2 Rabbit mAb                                                         |
| <b>Accession(Primary)</b>   | P09758                                                                   |
| <b>Host Species</b>         | Rabbit                                                                   |
| <b>Clonality</b>            | Monoclonal                                                               |
| <b>Isotype</b>              | IgG                                                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                        |
| <b>Calculated MW</b>        | Calculated MW: 36 kDa; Observed MW: 45-65 kDa                            |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:200-1:2000;<br>ICC/IF-1:50-1:100; FC-1:50-1:100 |

## Target information (P09758)

|                          |                                                                                                                                                                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | TACSTD2                                                                                                                                                                                                                                                                     |
| <b>About Gene Symbol</b> | TACSTD2 (Tumor Associated Calcium Signal Transducer 2)<br>is a protein kinase and signaling regulator involved in<br>adaptive immunity. It is an important research target in<br>cancer research and is widely studied to elucidate cellular<br>mechanisms and identify the |



Western blot analysis of TROP2 in HCT116, PC-3, MCF-7 lysates using TACSTD2 antibody.