

TAF1C Rabbit mAb


ABZ3184

Product Overview

Name	TAF1C Rabbit mAb
Accession(Primary)	Q15572
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 95 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q15572)

Gene Name	TAF1C
About Gene Symbol	TAF1C (TATA-Box Binding Protein Associated Factor, RNA Polymerase I Subunit C) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely s



Western blot analysis of TAF1C in K562 lysates
using TAF1C antibody.