

TAGAP Rabbit mAb

ABZ3185


Product Overview

Name	TAGAP Rabbit mAb
Accession(Primary)	Q8N103
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:100-1:200; IP-1:20-1:50

Target information (Q8N103)

Gene Name TAGAP

About Gene Symbol



Western blot analysis of TAGAP in rat Brain, C6, CHO-K1, HeLa lysates using TAGAP antibody.

