

Tau Rabbit mAb


ABZ3186

Product Overview


Name	Tau Rabbit mAb
Accession(Primary)	P10636
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 50-80 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:200-1:500

Target information (P10636)

Gene Name	MAPT
About Gene Symbol	



Western blot analysis of Tau in mouse brain, rat brain lysates using Tau antibody.



Western blot analysis of Tau in Human Hippocampus lysates using Tau antibody.