


## T Box Protein 3 Rabbit mAb

ABZ3188

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | T Box Protein 3 Rabbit mAb                 |
| <b>Accession(Primary)</b>   | O15119                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 79 kDa; Observed MW: 79 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                           |

### Target information (O15119)

**Gene Name** TBX3**About Gene Symbol**

Western blot analysis of Tbx3 in 3T3 lysates using  
T Box Protein 3 antibody.

