

TCTP Rabbit mAb


ABZ3191

Product Overview


Name	TCTP Rabbit mAb
Accession(Primary)	P13693
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:100-1:1000

Target information (P13693)

Gene Name	TPT1
About Gene Symbol	



Western blot analysis of TCTP in K562, rat Brain, C6, 3T3, Hela lysates using TCTP antibody.



Immunocytochemistry analysis of TCTP(green) in Hela using TCTP antibody, and DAPI(blue)